

DermaMax[®] RF

RF MICRONEEDLING SYSTEM

◆ Safe

◆ Exclusive

◆ Potent

◆ Compact

◆ Easy

◆ Modern




ESTETIQUE[®]
REINVENTING AESTHETICS

DermaMax[®] RF

RF MICRONEEDLING SYSTEM

Discover the latest in facial and body tightening and rejuvenation treatments with growing popularity in the aesthetic market: the new DermaMax RF Dermal Remodeling, an advanced high intensity fractional RF technology that effectively fights the signs of aging through a fast and safe treatment.

DermaMax RF is an aesthetic procedure that delivers high intensity fractional RF energy into the skin using sterilized microneedles. As they penetrate the tissues, these ultrafine needles produce microscopic thermal wounds at different depths inside the dermis to help create a new skin by activating its repair mechanisms. Furthermore, the RF energy stimulates natural collagen production, achieving a fabulous antiaging effect in minutes without any recovery time.



♦ Safe and effective for all skin types, from the lightest to the darkest.




MINIMAL
TO NO PAIN



NON INVASIVE



NO SIDE
EFFECTS



WITHOUT
RECOVERY TIME

DermaMax RF is a safe technique with minimum risk of discomfort that treats the skin from the inside out without damaging the epidermis, offering immediate results without side effects such as hyperpigmentation. The energy is delivered at different precise depths, right where it is needed, to provide a highly effective and safe treatment for all skin types and colors.

Multiply your profits with the best device in its category



- ♦ A powerful combination of microneedling and fractional RF energy for smooth, beautiful skin.

Microneedles

Its ultrafine needles penetrate up to 3.5 into the dermis and distribute the energy evenly in a focalized manner for noticeable and consistent results in every treatment session. The microscopic thermal wounds stimulate a repair response that helps create a new skin without damaging the tissues.

Fractional RF

The effect achieved by the microneedles, combined with the proven thermal properties of the high intensity fractional bipolar RF, provide an incredibly effective treatment by producing controlled heating that increases collagen production and reduces the signs of aging.

- ♦ Offer your clients a wide range of benefits in one single device

- Reduces fine lines and wrinkles
- Tightens the skin
- Lightens scars and stretchmarks
- Reduces spider veins
- Reduces enlarged pores
- Improves skin texture
- Achieves a more uniform skin tone

- ♦ Ultrafine needles minimize discomfort while offering an incredibly effective treatment

The science behind the new DermaMax RF

DermaMax RF's state-of-the-art applicator consists of sterilized microneedles that are applied to the treatment area, creating microscopic pathways within the dermis which remain open for 4 to 6 hours. These microneedles deliver high intensity fractional bipolar RF energy into the skin at 31 different precise depths to produce controlled heating without damaging the surrounding tissues.

The heat generated activates a natural body response that triggers the creation of new collagen and elastin fibers to restore the skin's youthful firmness, smoothness and elasticity.

- ♦ Achieve amazing results in just 3 sessions!

By the end of the treatment series the skin's appearance and texture significantly improve as it recovers its youthful firmness, smoothness, and elasticity. Fine lines and wrinkles fade, enlarged pores become smaller, the skin tone evens out, and scars and stretch marks become lighter, resulting in a new beautiful skin that feels and looks softer, brighter and healthier.

Rejuvenation
from the
inside out



♦ Works at a microscopic level to create a new younger, heathier looking skin.

While some microneedling systems can tear the skin, cause pain and get inconsistent results, the DermaMax RF technology works at different depths inside the skin to repair and rejuvenate it in a safe and consistent manner. Its ultrafine 0.2 mm needles minimize the risk of pain and do not cause bleeding, achieving visible results with minimum to no discomfort and without recovery time.

Its advanced design allows to adjust the settings at 31 different precise levels with a penetration of 0.5 to up to 3.5 mm in 0.1 mm increments, evenly distributing the heat for fast, exceptional results.



25 or 64 sterilized needles with gold plated tips penetrate up to 3.5 mm inside the skin for a deep, visible skin rejuvenation.

- Ergonomic design.
- Activation button incorporated on the handle for increased convenience.
- Easy to replace, disposable sterilized needle tips.
- Two needle tips for different treatment areas.
- Even pattern for consistent results.
- Biocompatible gold plating, suitable even for patients with metal allergies.
- The microneedles are insulated except for the tip, delivering the energy deep into the dermis without damaging the epidermis.

Applications	Technologies
<ul style="list-style-type: none">• Skin tightening• Reduction of fine lines and wrinkles• Revision of stretch marks and scars• Reduction of spider veins• Reduction of enlarged pores• Improvement in skin texture• More uniform skin tone	<ul style="list-style-type: none">• Microneedling• High intensity fractional bipolar RF

Estetique USA

2145 N Commerce Parkway
Weston, FL 33326 USA

Tel: 1 954-442-6595
Toll Free: 1 888 214 3856

Estetique Mexico

Avenida Homero, 1339 Piso 6
Colonia Polanco IV Sección
Miguel Hidalgo, CDMX CP
11550 Mexico

Tel: 52 55 4634 0496
52 55 5161 2483

E-mail: sales@estetique.com
ventas@estetique.com

Web: www.estetique.com

Más productos de Estetique

RejuveMax®
COLLAGEN LIGHT SYSTEM

LipoMax® RF
FACE | BODY

DermaMax®
HYDRODERMABRASION SYSTEM

LumiMax® PL
LASER SYSTEM

LumiMax®
LED THERAPY SYSTEM

LumiMax® CO2
CO2 LASER SYSTEM

- Two specially designed tips to treat any part of the body

25 needle tip

This tip accommodates perfectly to smaller areas including the eyes, the forehead and the neck, tightening the skin and effectively treating fine lines and wrinkles, enlarged pores and acne scars, among others.

64 needle tip

This size is ideal for body treatments such as skin tightening and stretch mark reduction in large areas such as the legs, the back, the abdomen, and the arms.

Features:

- Color 10.4 in touchscreen with user-friendly interface.
- Face and body pre-programmed treatments, and data saving button.
- Disposable, single use 0.2 mm high performance needles.
- 31 precise multilevel depth settings in 0.1 mm increments, from 0.5 to 3.5 mm.
- Two tips, a 25 needle one for small areas and a 64 needle one for larger areas.
- Ultrafine 0.2 mm needles for minimum to no pain.
- Highly biocompatible gold plated needles, safe for patients with metal allergies.
- Finger switch, foot switch and automatic mode.
- Safe for all skin tones including type VI, with minimal risk of post-inflammatory hyperpigmentation.

Technical Specifications

- | | |
|------------------|-------------------|
| • RF frequency | 1MHz |
| • Maximum energy | 35W |
| • Tip size | 5*5 mm and 8*8 mm |
| • Tip pins | 25 pins, 64 pins |
| • Needle size | 0.2 mm |
| • Depth | 0.5mm-3.5mm |
| • Screen | 10.4" touchscreen |
| • Dimensions | 50x50x110cm |
| • Net weight | 40kg |

